



## PRODUCT SPECIFICATION

catalog number:	<b>J10311</b>
product description:	Human IL-6 Antigen
lot number:	[lot specific]
source:	CHO
form:	Lyophilized
purity:	>95% by SDS-PAGE
buffer:	Lyophilized from PBS
preservative:	None
molecular weight:	21-23 kDa by SDS-PAGE
endotoxin:	< 0.2 EU/μg by LAL method
reconstitution:	Reconstitute in ddH <sub>2</sub> O or PBS at 100 μg/ml
storage:	Lyophilized -20°C (up to 6 months from date of receipt); After reconstitution: 2-8°C (one month); -20°C (up to 2 months)
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.