



PRODUCT SPECIFICATION

catalog number:	J10312
product description:	Human IL-6 Antigen, Recombinant
lot number:	[lot specific]
source:	<i>E. Coli</i>
form:	Liquid; aseptically filtered
purification method:	Ion exchange chromatography
purity:	≥97%
buffer:	10 mM sodium phosphate, 150mM NaCl, 25% glycerol, 1mM AEBSE, pH 7.4 ± 0.2
preservative:	None
molecular weight:	20.9 kDa
concentration:	[lot specific] mg/mL (by ECLIA)
bioburden:	≤ 10 CFU/mL
manufacture date:	[lot specific]
expiration date:	[lot specific]
storage:	-20°C or colder
hazard/biohazard:	There is no known hazard or biohazard associated with this product.