



PRODUCT SPECIFICATION

catalog number:	J10315566
product description:	Human IL-6 Protein, Recombinant
lot number:	[lot specific]
source:	<i>E. coli</i> , Pro29-Met212
form:	Liquid; 0.2um filtered
concentration:	[lot specific]
purity:	>97% by SDS-PAGE under reducing conditions
buffer:	PBS + 1M NaCl, pH 7.4 with 50ug BSA per 1ug as carrier protein
preservative:	None
molecular weight:	~20.3 kDa
storage:	Upon receipt store -20°C to -70°C (up to 12 months)
hazard/biohazard:	There is no known hazard or biohazard associated with this product.