



PRODUCT SPECIFICATION

CATALOG NUMBER: J10700

PRODUCT DESCRIPTION:	Recombinant Human Pro-Brain natriuretic peptide (rNT-ProBNP)
STRUCTURE:	aa 1-76 of N-terminal Pro-BNP + His Tag
SOURCE:	<i>E. coli</i>
BUFFER:	50 mM potassium phosphate buffer, 250mM NaCl, pH 6.7
PRESERVATIVE:	none
PURITY:	97.8% by SDS-PAGE; density scanning
MOLECULAR WEIGHT:	8.5 kD by SDS-PAGE
STORAGE:	-20°C. Avoided repeated freeze thaw cycles.
HAZARD/BIOHAZARD:	There is no known hazard or biohazard associated with this product.