



## PRODUCT SPECIFICATION

CATALOG NUMBER: J10710

PRODUCT DESCRIPTION:	Recombinant Human Pro-Brain natriuretic peptide (rProBNP)
STRUCTURE:	aa 1-108 of Pro-BNP + His Tag
SOURCE:	<i>E. coli</i>
BUFFER:	20 mM Tris HCl, 500 mM NaCl, pH 8.0 with protease inhibitors
PRESERVATIVE:	none
PURITY:	90% by SDS-PAGE; density scanning
MOLECULAR WEIGHT:	13.7 kD by SDS-PAGE
STORAGE:	-20°C. Avoided repeated freeze thaw cycles.
HAZARD/BIOHAZARD:	There is no known hazard or biohazard associated with this product.