



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J11010**

PRODUCT DESCRIPTION:	TSH Antigen, standard grade
SOURCE:	Human Pituitary Glands
FORM:	Lyophilized powder, lyophilized from 50mM ammonium bicarbonate
PURITY:	>50% by SDS-PAGE
ACTIVITY:	4.0 IU/mg
CONTAMINANTS:	LH: 0.7021% FSH: 0.00689%
STORAGE TEMPERATURE:	2-8°C short term -20°C long term
HAZARD/BIOHAZARD:	Donors tested negative for anti-HIV 1 & 2, anti-HCV, Syphilis, HBcAg, and HBsAg by currently approved FDA methods. There is no known hazard associated with this product.