



PRODUCT SPECIFICATION

CATALOG NUMBER: J11030

PRODUCT DESCRIPTION:	TSH Antigen, Recombinant
SOURCE:	Chinese Hamster Ovary Cells.
FORM:	Lyophilized
PURITY:	≥95% by SDS-PAGE
FILTERED:	0.2 µm filter
ACTIVITY:	9 IU/mg
STORAGE TEMPERATURE:	-20°C long term. After reconstitution, store material at 2-8°C for one month or -20°C for three months in manual defrost freezer. Avoid repeated freeze-thaw cycles.
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.