



PRODUCT SPECIFICATION

CATALOG NUMBER: J11120

PRODUCT DESCRIPTION:	TSH Antigen, Beta Subunit
SOURCE:	Human Pituitary Glands
FORM:	Lyophilized
PURITY:	>98% by SDS-PAGE
CONTAMINANTS:	TSH 1.247% FSH 0.002% LH 0.013%
STORAGE TEMPERATURE:	2-8°C short term, -20°C long term
HAZARD/BIOHAZARD:	Serum from donors has been tested and found negative for anti-HCV, HBsAg, HBcAg, anti-HIV 1 and 2, and Syphilis by currently approved FDA methods. There is no known hazard associated with this product.