



## PRODUCT SPECIFICATION

catalog number:	<b>J12020</b>
product description:	Human Luteinizing Hormone (LH) Antigen, High Purity
lot number:	[lot specific]
source:	Human tissue
form:	Lyophilized powder
purity:	>95%
buffer:	Lyophilized from 10mM NH <sub>4</sub> HCO <sub>3</sub>
activity:	[lot specific] IU/mg by WHO 1 <sup>st</sup> IRP68/40
manufacture date:	[lot specific]
expiration date:	3 years from date of manufacture
reconstitution:	Reconstitute in a phosphate buffered saline with 0.5% BSA and store in small aliquots at -20°C ; avoid multiple freeze thaw cycles.
storage:	2-8°C as lyophilized
hazard/biohazard:	Donor material has been tested and found negative for HIV, HBsAg and HCV by FDA approved methods. However, this material should still be handled with caution as potentially hazardous.