



PRODUCT SPECIFICATION

CATALOG NUMBER: J12050

PRODUCT DESCRIPTION:	LH Antigen, Affinity Pure
SOURCE:	Human Pituitary Glands
FORM:	Lyophilized powder, lyophilized from 50mM ammonium bicarbonate, no preservative added.
PURITY:	98% by SDS-PAGE
REFERENCE ASSAY:	Abbott IMX
REFERENCE PREP:	2nd IS 80/552
CONTAMINANTS:	FSH TSH HGH HPRL
STORAGE TEMPERATURE:	2-8°C lyophilized. Once reconstituted, dispense in small aliquots and store at -20°C.
BIOHAZARD:	Donor tested negative for HIV 1&2, HCV, Syphilis, HBcAg, and HBsAg by FDA approved methods.