



PRODUCT SPECIFICATION

CATALOG NUMBER: J15120

PRODUCT DESCRIPTION:	Prolactin Antigen, Standard grade
SOURCE:	Human pituitary glands
FORM:	Lyophilized powder; lyophilized from 50 mM ammonium bicarbonate
PURITY:	>50% by SDS-PAGE
SPECIFIC ACTIVITY:	14.3 IU/mg
REFERENCE ASSAY:	DPC Immulite 1000
REFERENCE PREPARATION:	THIRD IS (84/500), WHO
CONTAMINANTS:	LH 0.021% HG 0.32%
STORAGE TEMPERATURE:	2-8°C short term -20°C long term
HAZARD/BIOHAZARD:	Serum from donors has been tested and found to be negative for HIV-1 and 2, Hepatitis B Core antigen, Hepatitis B Surface antigen, and Hepatitis C virus by current approved FDA methods. There is no known hazard associated with this product.