



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J15130**

PRODUCT DESCRIPTION:	Prolactin, Recombinant Human Antigen
SOURCE:	<i>Escherichia Coli</i>
FORM:	Lyophilized powder; lyophilized from a 1 mg/ml solution in 10 mM sodium phosphate buffer, 50 mM NaCl, pH 8.0. Reconstitute in sterile deionized water at 0.1 to 1 mg/ml.
PURITY:	>95% by SDS-PAGE
STORAGE TEMPERATURE:	2-8°C short term, -20°C long term. After reconstitution, store at -20°C. Avoid repeat freeze thaw cycles.
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.