



PRODUCT SPECIFICATION

CATALOG NUMBER: J15140

PRODUCT DESCRIPTION:	Human Prolactin Antigen, Recombinant
SOURCE:	<i>E. coli</i>
FORM:	Liquid
BUFFER:	20 mM Tris-HCl, 10 mM EDTA, pH 8
PURITY:	90% (by SDS-PAGE)
MOLECULAR WEIGHT:	24 kDa
STORAGE:	2-8°C short term, -20°C long term; avoid repeated freeze/thaw cycles
HAZARD/BIOHAZARD:	There is no known hazard or biohazard associated with this product.