



## PRODUCT SPECIFICATION

catalog number:	<b>J15150</b>
product description:	Prolactin Antigen, High Pure
lot number:	[lot specific]
source:	Human
form:	Lyophilized
purity:	>95% Iodination grade
activity:	[lot specific]
cross-reactivity:	<ul style="list-style-type: none"><li>▪ hGH                    0.05%</li><li>▪ hLH                    0.100%</li><li>▪ hTSH                  0.001%</li></ul>
reconstitution:	Reconstitute in PBS with 0.5% BSA.
storage:	2-8°C lyophilized; -20°C after reconstitution
hazard/biohazard:	Human source has been tested and found negative for HIV, HBsAg and HCV. However, this material should still be handled with caution as potentially hazardous.