



**PRODUCT SPECIFICATION**

**CATALOG NUMBER: J16000**

PRODUCT DESCRIPTION:	HGH Antigen, Affinity Purified
SOURCE:	Human Pituitary Glands
FORM:	Lyophilized powder
PURITY:	>98%; as determined by SDS-PAGE
SPECIFIC ACTIVITY:	2.25 IU/mg
REFERENCE PREPARATION:	First IS (80/505), WHO
CONTAMINANTS:	hTSH < 0.003% hLH < 0.003% hFSH < detectable limit hPRL < 0.1
STORAGE TEMPERATURE:	2-8°C short term, -20 <sup>0</sup> C long term
VIRAL TESTS:	Product tested negative for HBsAg, HIV 1&2, and HCV by FDA approved test methods.