



PRODUCT SPECIFICATION

catalog number:	J18050
product description:	Human Follicle Stimulating Hormone (FSH) Antigen, Recombinant
lot number:	[lot specific]
source:	Chinese hamster ovary cells
form:	Lyophilized powder
purity:	>95% by SDS-PAGE
buffer:	Lyophilized from a 0.2 μ filtered solution in PBS with 5% Trehalose, pH 7.2.
activity:	[lot specific] IU/mg (by Abbott Architect)
reconstitution:	Reconstitute with a sterile PBS to 100 μ g/mL.
storage:	2-8°C short term (one month); -20°C long term avoiding multiple freeze thaw cycles.