



PRODUCT SPECIFICATION

CATALOG NUMBER: J19300

PRODUCT DESCRIPTION:	Thyroxine Binding Globulin Antigen
SOURCE:	Human serum
FORM:	Lyophilized from 0.02M NH_4HCO_3 . May contain traces of buffer salts.
PURITY:	>98% by SDS-PAGE (one band at 55kD)
PRESERVATIVE:	None
RECONSTITUTION:	PBS, pH>7.0 containing 0.15M NaCl
STORAGE TEMPERATURE:	2-8°C
HAZARD/BIOHAZARD:	Donor serum tested and found negative for anti-HIV 1/2, anti-HCV and HBsAg antigen. There is no known hazard associated with this product.