



## PRODUCT SPECIFICATION

catalog number:	J29350601
product description:	Calprotectin Dimeric Recombinant Protein (S100A8 & S100A9)
lot number:	[lot specific]
source:	E.coli
form:	Liquid
buffer:	Tris buffer, pH 8.0
preservative:	0.09% Sodium Azide
concentration:	[lot specific] mg/ml
expiration date:	[lot specific]
storage:	-20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known hazard or biohazard associated with this product.