



## PRODUCT SPECIFICATION

catalog number:	J29352601
product description:	Calprotectin Tetrameric Recombinant Protein
lot number:	[lot specific]
source:	E.coli
form:	Liquid
buffer:	Tris buffer, pH 8.0
preservative:	0.09% Sodium Azide
concentration:	[lot specific] mg/ml
specificity:	Recombinant tetrameric calprotectin (2 subunits S100A8 and 2 subunits S100A9) modified to maintain tetrameric conformation
expiration date:	[lot specific]
storage:	-20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known hazard or biohazard associated with this product.