



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J31520**

PRODUCT DESCRIPTION:	Recombinant human Creatine Kinase MM Isoenzyme (CKMM) Antigen.
SOURCE:	<u>E. coli</u>
FORM:	Liquid
PURITY:	>90% by SDS-PAGE. 45kD protein.
BUFFER:	50 mM Tris-HCl buffer, 5 mM beta mercaptoethanol, pH 8.0, 1 mM EDTA, 5 mM disodium succinate
STORAGE TEMPERATURE:	-20°C or lower.