



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J33000**

PRODUCT DESCRIPTION:	Recombinant human ventricular myosin light chain-1
SOURCE:	E.coli
PURITY:	90% by SDS-PAGE (density scanning)
BUFFER:	50 mM Tris-HCL, 1 mM DTT, 200mM NaCl, pH 8.0
MOLECULAR WEIGHT:	28 kD by SDS-PAGE
STORAGE TEMPERATURE:	-20°C long term
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.