



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J34010**

PRODUCT DESCRIPTION:	Troponin I Antigen, (Human Skeletal Muscle Troponin I)
SOURCE:	Human skeletal muscle
PURITY:	≥98% by SDS-PAGE: ~21-22 kDa band with minor fragmentation bands
BUFFER:	50 mM Tris-HCl, 100 mM NaCl, 8 M Urea, 15 mM β-mercaptoethanol, 2 mM EDTA, pH 7.4
STORAGE TEMPERATURE:	-20°C
HAZARD/BIOHAZARD:	Sample from tissue donor was tested and found negative for HBsAg, HIV-1 and HIV-2 antibodies, HCV, and syphilis. There is no known hazard associated with this product.