



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J34120**

PRODUCT DESCRIPTION:	Recombinant Human Cardiac Troponin I
SOURCE:	Expressed in <i>E. coli</i>
FORM:	Liquid
PURITY:	≥98% by SDS PAGE. Major 24 kDa band.
BUFFER:	6M Urea, 50 mM Tris, pH 8.0
STORAGE TEMPERATURE:	-20°C long term (avoid repeated freeze/thaw cycles)
HAZARD/BIOHAZARD:	There is no known hazard or biohazard associated with this product.