



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J34121**

PRODUCT DESCRIPTION:	Recombinant Troponin I Antigen, High Pure
SOURCE:	<i>E. Coli</i>
FORM:	Liquid
PURITY:	≥98% by SDS-PAGE <sup>12%</sup> Major band at 26.1 kD
PURIFICATION:	Metal affinity chromatography
BUFFER:	20mM Tris, 150mM NaCl, 8M Urea, pH 8.0±0.2
PRESERVATIVE:	None
BIOBURDEN:	< 10 CFU/ml
STORAGE TEMPERATURE:	-20°C (Avoid repeated freeze/thaw cycles)
HAZARD/BIOHAZARD:	There is no known hazard or biohazard associated with this product.