



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J34122**

PRODUCT DESCRIPTION:	Stabilized Recombinant Human Cardiac Troponin I Antigen
SOURCE:	Expressed in E. coli.
FORM:	Liquid
BUFFER:	500 mM NaCl, 5 mM sodium phosphate, 2 mM CaCl <sub>2</sub> , 0.5 mM EDTA, (pH 7.0) containing stabilizer
PRESERVATIVE:	0.05% NaN <sub>3</sub>
STORAGE TEMPERATURE:	-20 to -70°C (Avoid repeated freeze/thaw cycles)
COMMENTS:	Prepare stock dilutions in delipidated human serum containing 1 mM CaCl <sub>2</sub> . Freeze as aliquots and store at -70°C.
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.