



PRODUCT SPECIFICATION

CATALOG NUMBER: J34300

PRODUCT DESCRIPTION:	Fatty Acid Binding Protein Antigen
SOURCE:	Human Heart Tissue
FORM:	Lyophilized powder
PURITY:	>95% by SDS-PAGE
BUFFER:	10 mM Tris HCl, pH 8.0
RECONSTITUTION:	Recommended to reconstitute with distilled water
STORAGE TEMPERATURE:	2-8°C short term, -20°C long term
HAZARD/BIOHAZARD:	Blood sample from tissue donor teated and found negative for HBsAg, HIV-1/2 antibodies, HCV antibody, and syphilis. There is no known hazard associated with this product.