



PRODUCT SPECIFICATION

CATALOG NUMBER: J34510

PRODUCT DESCRIPTION:	Recombinant Cardiac Troponin T (rTnT)
SOURCE:	<i>E. coli</i>
FORM:	Liquid
PURITY:	>90%
BUFFER:	50 mM Tris-HCl, 500 mM NaCl, 2 mM EDTA, 15 mM beta-mercaptoethanol, pH 7.5
PRESERVATIVES:	None
STORAGE TEMPERATURE:	-20°C
HAZARD/BIOHAZARD:	There is no known hazard or biohazard associated with this product.