



## PRODUCT SPECIFICATION

CATALOG NUMBER: J35000

|                      |   |
|----------------------|---|
| PRODUCT DESCRIPTION: | Purified Human Troponin T Skeletal Muscle Antigen   |
| SOURCE:              | Human skeletal muscle   |
| PURITY:              | > 95% by SDS-PAGE   |
| FORM:                | Lyophilized powder, from 0.01 M HCl   |
| RECONSTITUTION:      | Reconstitute with distilled water   |
| PRESERVATIVE:        | None  |
| STORAGE:             | -20°C short term; -70°C long term   |
| NOTE:                | Blood sample from tissue donor was tested and found negative for HBsAg, HIV 1&2 antibody, and HCV antibody. There is no known biohazard associated with this product. |