



PRODUCT SPECIFICATION

CATALOG NUMBER: J36000

PRODUCT DESCRIPTION:	Human Cardiac Troponin C Recombinant Antigen
FORM:	Liquid
PURITY:	95%, by SDS PAGE
BUFFER:	150 mM NaCl, 10 mM sodium phosphate, pH 7.0
PRESERVATIVE:	0.05 NaN ₃
STORAGE TEMPERATURE:	-20°C (Avoid repeated freeze/thaw cycles)
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.