



PRODUCT SPECIFICATION

CATALOG NUMBER: J36010

PRODUCT DESCRIPTION:	Recombinant Cardiac Troponin C
SOURCE:	Expressed in <u>E coli</u> .
FORM:	Liquid
BUFFER:	50 mM Tris-HCl, 0.5 mM NaCl, 2mM EDTA, 5 mM Beta-Mercaptoethanol, pH 7.5
STORAGE TEMPERATURE:	-20°C
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.