



PRODUCT SPECIFICATION

CATALOG NUMBER: J36050

PRODUCT DESCRIPTION:	Galectin-3 Recombinant Antigen (with HisTag)
SOURCE:	<i>E. Coli</i>
FORM:	Liquid
APPEARANCE:	Clear to colorless
PURIFICATION:	Metal affinity chromatography
PURITY:	>95% by SDS-PAGE ^{12%} Major band at 31.0 kD
BUFFER:	Asceptically filtered liquid in 20mM Tris-HCL, 0.1 mM DTT, 10% glycerol, pH 8.0 ± 0.2
BIOBURDEN:	<10 CFU/ml
STORAGE:	-20°C
HAZARD/BIOHAZARD:	There is no known hazard or biohazard associated with this product.