



PRODUCT SPECIFICATION

catalog number:	J36330
product description:	HbA1c Antigen, High Pure
lot number:	[lot specific]
source:	Human erythrocytes
form:	Clear, red frozen solution; 0.2 µm filtered
purity:	>96% by HPLC
buffer:	Buffer formulation, pH 8.0
preservative:	None
concentration:	[lot specific] mg/mL (determined by OD _{280nm} at 415nM, E=7)
bioburden:	< 100 CFU/mL
manufacture date:	[lot specific]
storage:	-20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	Donor material has been tested and found negative for HBsAg, anti-HIV 1&2, and anti-HCV by FDA approved methods. However, this material should still be handled with caution as potentially hazardous.