



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J37020**

PRODUCT DESCRIPTION:	Glycogen Phosphorylase BB, Human heart
SOURCE:	Human heart
FORM:	Liquid
PURITY:	>95% (SDS-PAGE)
BUFFER:	1 mM Beta-glycerophosphate, 1mM EDTA, 15 mM ME, 0.5 M NaCl, pH 7.8, (50% glycerol has been added)
STORAGE TEMPERATURE:	2-8°C short term, -20°C long term
BIOHAZARD:	From tissue of individuals that have been shown by certified tests to be negative for HBsAg, HCV and for HIV 1&2 antibodies.