



PRODUCT SPECIFICATION

CATALOG NUMBER: J38030

PRODUCT DESCRIPTION:	Human D-Dimer Antigen, Part Pure
SOURCE:	Human plasma
FORM:	Liquid
BUFFER:	50mM phosphate buffer, 0.15M sodium chloride, pH 7.8 ± 0.2, 20% glycerol
PRESERVATIVE:	0.09% NaN ₃
CONCENTRATION:	Reference (by ELFA method on VIDAS)
TOTAL PROTEIN:	Reference (by Lowry & Bradford method)
STORAGE:	2-8°C
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precaution. There is no known biohazard associated with this product