



PRODUCT SPECIFICATION

CATALOG NUMBER: J41420

PRODUCT DESCRIPTION:	Mumps Antigen
FORM:	Liquid
STRAIN:	Enders (ATCC VR-106)
PURIFICATION:	Vero cells are infected with the Enders strain of Mumps virus. Antigen is prepared from the culture supernatant of the infected cell cultures. The antigen preparation is partially purified to enrich for viral antigen and reduce host cell protein.
BUFFER:	60 mM glycine, NaCl buffer, pH 9.5
STORAGE TEMPERATURE:	-70°C
APPLICATION:	Antigen is intended for use in solid phase immunoassays, giving high sensitivity and low background. For <u>in vitro</u> use only.
HAZARD/BIOHAZARD:	This material contains Sodium Azide as preservative, take appropriate precautions. Although the antigen has been tested on Vero cells and produces no CPE, the antigen should be handled as if potentially infectious.