



PRODUCT SPECIFICATION

CATALOG NUMBER: J41610

PRODUCT DESCRIPTION:	Varicella Zoster antigen
STRAIN:	Wild strain isolate from skin lesion of a child with Varicella.
FORM:	Liquid frozen at -70°C
PURIFICATION:	Antigen is prepared by alkaline extraction of human fibroblast cells infected with Varicella Zoster virus. It is partially purified to reduce host cell components.
BUFFER:	60 mM glycine buffer pH 9.6
PRESERVATIVE:	0.1% NaN_3
STORAGE TEMPERATURE:	-70°C
APPLICATION:	Antigen performs well in solid phase immunoassays, giving high sensitivity and a low background. For in vitro use only.