



PRODUCT SPECIFICATION

CATALOG NUMBER: J41630

PRODUCT DESCRIPTION:	Recombinant VZV ORF9 protein
SOURCE:	E.coli derived recombinant. The protein contains the VZV ORF9 immunodominant regions.
FORM:	Liquid
PURIFICATION METHOD:	Sepharose-derived purification
PURITY:	>95% by 10 % PAGE (coumassie staining)
BUFFER:	25mM Tris-HCl, 1mM EDTA, 50% glycerol
STORAGE TEMPERATURE:	Shipped at ambient temperature. Upon arrival store at -20 °C.
STABILITY:	Five years frozen; one month in solution at room temperature
SPECIFICITY:	Immunoreactive with sera of VZV-infected individuals
APPLICATIONS:	ELISA, Western blot