



PRODUCT SPECIFICATION

catalog number:	J43050
product description:	Cytomegalovirus (CMV) Pp150 Recombinant Antigen, aa 1011-1048
lot number:	[lot specific]
source:	<i>E. coli</i>
form:	Liquid
purification method:	GS-4B sepharose-affinity purification
purity:	>95% pure by SDS-PAGE (Coomassie staining)
buffer:	25 mM Tris-HCl, 1 mM EDTA, 50% glycerol, pH 8.0
isotype:	IgG ₁
concentration:	[lot specific] mg/mL
specificity:	Immunoreactive with sera of CMV-infected individuals.
application/notes:	This antibody is stable at room temperature in solution for one month or for up to five years if kept frozen.
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.