



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J43070**

PRODUCT DESCRIPTION:	Cytomegalovirus (CMV) Pp38 recombinant antigen aa 117-373
SOURCE:	<u>E. coli</u>
PURIFICATION METHOD:	GS-4B Sepharose-affinity purification
PURITY:	>95% pure by SDS-PAGE (Coomassie staining)
BUFFER:	25 mM Tris-HCl, 1mM EDTA, 50% glycerol, pH 7.2
STORAGE TEMPERATURE:	2-8 °C short term, -20°C long term
STABILITY:	Five years frozen
SPECIFICITY:	Immunoreactive with sera of CMV-infected individuals.
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.