



PRODUCT SPECIFICATION

CATALOG NUMBER: J43500

PRODUCT DESCRIPTION:	EBV Antigen Cell Extract
SOURCE:	P3H3 Human Burkitt's Lymphoma cells
FORM:	Liquid
BUFFER:	Tris containing NP40, sodium deoxycholate and PMSF
PRESERVATIVE:	None
STORAGE TEMPERATURE:	-70°C
PREPARATION:	Chronically infected cells are harvested by centrifugation and disrupted by sonication and cell debris is removed by centrifugation.
BIOHAZARD:	All materials should be handled as if potentially infectious. Generally accepted laboratory practices for infectious materials should be employed when handling this product.