



## PRODUCT SPECIFICATION

catalog number:	<b>J43510149</b>
product description:	EBV (Epstein Barr Virus) Viral Capsid gp125 Antigen
lot number:	[lot specific]
form:	Liquid
preparation:	Chronically infected cells are harvested by centrifugation and disrupted by sonication and cell debris is removed by centrifugation. This lysate is purified on an immunoaffinity column using a monoclonal antibody against gp125. .
buffer:	20mM Tris-HCl, 3M MgCl <sub>2</sub> , pH 7.4
concentration:	[lot specific] µg/mL (by Biorad dye binding assay)
storage:	-70°C to -100°C avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product. However, this material should still be handled with caution as potentially hazardous.