



PRODUCT SPECIFICATION

CATALOG NUMBER: J43520

PRODUCT DESCRIPTION:	EBV VCA gp125 Antigen. Epstein Barr Virus, Viral Capsid Antigen, gp 125
FORM:	Liquid
PURIFICATION:	Antigen was affinity purified using a monoclonal antibody specific for VCA gp125 protein.
BUFFER:	0.02M Tris containing 3M Magnesium Chloride, pH 6.4-7.2
PRESERVATIVE:	None
STORAGE TEMPERATURE:	-70°C
INACTIVATION:	Sodium deoxycholate and NP-40 with sonication
BIOHAZARD:	All materials should be handled as if potentially infectious. Generally accepted laboratory practices appropriate for infectious materials should be employed when handling this product.