



PRODUCT SPECIFICATION

catalog number:	J43540
product description:	EBV Mosaic EBNA1 Antigen, Recombinant
lot number:	[lot specific]
source:	<i>E. coli</i> derived; contains the HHV-4 EBNA region
form:	Liquid
purification method:	Sepharose-derived purification
purity:	>95% by SDS-PAGE (coomassie staining)
buffer:	50mM Tris-HCl, 10mM glutation, 60mM NaCl, 0.5% sarcosyl, pH 8.0
concentration:	[lot specific] mg/mL
specificity:	This protein is immunoreactive with sera of EBV-infected individuals.
application/notes:	Stability has been shown for up to one month at room temperature in solution or for up to five years if kept frozen.
storage:	Shipped at 2-8°C; upon arrival store at -20°C avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.