



PRODUCT SPECIFICATION

CATALOG NUMBER: J43550

PRODUCT DESCRIPTION:	EBV p18 Recombinant Antigen
SOURCE:	<u>E. coli</u>
PURIFICATION METHOD:	Sepharose-derived purification
PURITY:	95% pure by SDS-PAGE (Coomassie staining)
BUFFER	25 mM Tris-HCl, 1.5 M urea, 50% glycerol
STORAGE TEMPERATURE:	-20°C
SPECIFICITY:	Immunoreactive with sera of EBV-infected individuals.
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.