



PRODUCT SPECIFICATION

| | |
|----------------------|--|
| catalog number: | J43610 |
| product description: | Influenza A (H1N1) Antigen |
| lot number: | [lot specific] |
| source: | Cultured in embryonated chicken eggs; H1N1 subtype |
| form: | Liquid |
| buffer: | PBS |
| concentration: | [lot specific] mg/mL |
| titer: | [lot specific] HA units/mL by Hemagglutination Endpoint Assay |
| storage: | -70°C to -100°C avoiding multiple freeze thaw cycles. |
| inactivation: | Virus preparation is inactivated by Gamma radiation. Effectiveness is determined when no detection of growth in two blind passages under the optimal culture conditions. |
| hazard/biohazard: | This material has been inactivated and should contain no infectious material. However, this material should still be handled with caution as potentially hazardous. |