



PRODUCT SPECIFICATION

catalog number:	J43617
product description:	Influenza A Antigen, Recombinant
lot number:	[lot specific]
source:	E. Coli
form:	Clear Liquid
purity:	>90%
buffer:	Tris Buffer with additives and salts, pH 8.0
preservative:	0.09% NaN ₃
concentration:	[lot specific] mg/ml
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles
expiration date:	[lot specific]
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.