



PRODUCT SPECIFICATION

CATALOG NUMBER: J43620

PRODUCT DESCRIPTION:	Influenza B Antigen
STRAIN:	Hong Kong 5/72
BUFFER:	EMEM
TITER:	40,960 HA units/ml by Hemagglutination Endpoint Assay
STORAGE TEMPERATURE:	<-70°C; avoid repeated freezing and thawing
PREPARATION:	Clarified and concentrated cell culture supernatant.
INACTIVATION:	Virus preparation is inactivated by Gamma radiation.
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.