



#### PRODUCT SPECIFICATION

catalog number:	<b>J44015126</b>
product description:	HBsAG Antigen, Subtype adr, recombinant
lot number:	J0198
source:	Yeast
form:	Liquid
purity:	>98% (by SDS-PAGE)
buffer:	10mM PO <sub>4</sub> , 140mM Sodium Chloride, pH 6.4
preservative:	0.1% NaN <sub>3</sub>
concentration:	1.27 mg/ml
storage:	2-8°C; DO NOT FREEZE
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.