



PRODUCT SPECIFICATION

catalog number:	J44200
product description:	Hepatitis B Virus E (HBeAg), Recombinant Antigen
lot number:	[lot specific]
source:	<i>E. coli</i> ; contains the immunodominant region
form:	Liquid
purification method:	GS-4B sepharose-affinity purification
purity:	>95% pure by SDS-PAGE (Coomassie staining)
buffer:	25mM Tris-HCl, 1.5M urea, 50% glycerol, pH 8.0
concentration:	[lot specific] mg/mL
specificity:	Immunoreactive with sera of HBV-infected individuals.
application/notes:	This protein is stable at room temperature in solution for one month or for up to five years if kept frozen.
storage:	2-8°C short term (one month); upon arrival store at -20°C long term (five years), avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.