



PRODUCT SPECIFICATION

CATALOG NUMBER: J44300

| | |
|----------------------|--|
| PRODUCT DESCRIPTION: | Hepatitis A Antigen |
| FORMULATION: | Prepared from the supernatant of infected Wannenstapel and dialysed against PBS. |
| BUFFER: | Phosphate buffered saline |
| STORAGE TEMPERATURE: | -20°C |
| INACTIVATION: | Inactivated by incubating with formaldehyde for 11 days. |
| HAZARD/BIOHAZARD: | This material contains sodium azide as a preservative, take appropriate precautions. The reagent has been inactivated and should contain no infectious material. There is no known biohazard associated with this product. |